

*Insertion via mini-port (Covidien Japan); insertion using a shaft guide (lower image)

MANUFACTURER

Hirata Precisions Co., Ltd.

2-10-17 Kunugiyama Kamagaya-city, Chiba, JAPAN.

TEL 81-47-386-2101

E-mail: hiratapr@bolero.plala.or.jp

Endo Relief Small-Diameter Forceps 12B2X00014000028 Vender

SALES

Hope Denshi Co., Ltd.

2-10-17 Kunugiyama Kamagaya-city, TEL 81-47-389-5755 FAX 81-47-386-8645 http://www.hope-denshi.co.jp E-Mail: info@hope-denshi.co.jp

For the surgical era of reduced port and single port access

Endo Relief

Advanced Endo Relief forceps for unrivalled operating comfort

Needle-like forceps

Endo Relief

Endo Relief forceps provide capability far beyond the realm of conventional forceps. The forceps tips are the same design as the 5-mm tips used in conventional laparoscopic surgery and can therefore be used for cholecystectomy, gastrointestinal procedures, as well as other forms of laparoscopic surgery. With comfort like this, as you monitor your actions onscreen, it's easy to forget that the shaft has a diameter of only 2.4 mm!

Fine craftsmanship enables the handle of the 2.4-mm forceps to be detached and reattached.

The small-diameter trocar mini-port leaves a virtually invisible wound.



Ergonomic design

Smooth curvature extends from the 5-mm tips to the 2.4-mm diameter shaft.

Sensitive Operability The tips can be rotated in the conventional manner through

handle-side control. With such smooth operation, it's hard to believe the small diameter of these forceps.

Rigidity and Durability

A broad insertion section with a titanium alloy shaft provides high rigidity and exceptional durability.





Safety and Reuse

To deliver safety through a familiar feel and enable repeated sterilization, the tips have been designed with the same size and form as conventional 5-mm forceps.

